


Product Name: Recombinant Human XBP1 Protein
Catalog #: PEH3534



Summary

Name	XBP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17861
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQRLTHLSPE EKALRRKLKN RVAAQTARDR KKARMSELEQ QVVDLEENQ KLLLENQLLR EKTHGLVVEN QELRQRLGMD ALVAEEEAEG KGNEVRPVAG

image



Background

Alternative Names	X-box-binding protein 1 (XBP-1) (Tax-responsive element-binding protein 5) (TREB-5) [Cleaved into: X-box-binding protein 1, cytoplasmic form; X-box-binding protein 1, luminal form]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use