


Product Name: Recombinant Human NCF2 Protein
Catalog #: PEH3535



Summary

Name	NCF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P19878
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSPGQKQKEE PKEVKLSVPM PYTLKVHYKY TVVMKTQPGL PYSQVRDMVS KKLELRLEHT KLSYRPRDSN ELVPLSEDSM KDAWGQVKNY CLTLWCENTV GDQGFPDEPK ESEKADANNQ TTEPQLKKGS QVEALFSYEA TQPEDLEFQE GDIIVLVSKV NEEWLEGECK GKVGIKPKVF VEDCATTDLE

image



Background

Alternative Names	Neutrophil cytosol factor 2 (NCF-2) (67 kDa neutrophil oxidase factor) (NADPH oxidase activator 2) (Neutrophil NADPH oxidase factor 2) (p67-phox)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NCF2 Protein
Catalog #: PEH3535

