


Product Name: Recombinant Human TRH Protein
Catalog #: PEH3537



Summary

Name	TRH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20396
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLFLRENIQ RLQGDQGEHS ASQIFQSDWL SKRQHPGKRE EEEEEGVEEE EEEEGGAVGP HKRQHPGRRE DEASWSVDVT QHKRQHPGRR

image



Background

Alternative Names	Pro-thyrotropin-releasing hormone (Pro-TRH) (Prothyroliberin) [Cleaved into: Thyrotropin-releasing hormone (TRH) (Protirelin) (TSH-releasing factor) (Thyroliberin) (Thyrotropin-releasing factor) (TRF)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use