


Product Name: Recombinant Human ATP2B4 Protein
Catalog #: PEH3542



Summary

Name	AT2B4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P23634
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QCRIEQEQKF SIIRNGQLIQ LPVAEIVVGD IAQVKYGDLL PADGILIQGN DLKIDESSLT GESDHVKKSL DKDPMLLSGT HVMEGSGRMV VTAVGVNSQT GIILTLLGVN EDDEGEKKKK GKKQGVPENR NKAKTQDGVA LEIQPLNSQE GIDNEEKDKK AVKVPKKEKS

image



Background

Alternative Names	Plasma membrane calcium-transporting ATPase 4 (PMCA4) (EC 7.2.2.10) (Matrix-remodeling-associated protein 1) (Plasma membrane calcium ATPase isoform 4) (Plasma membrane calcium pump isoform 4)
--------------------------	--

Note

Product Name: Recombinant Human ATP2B4 Protein
Catalog #: PEH3542



For Research Use Only , Not for Diagnostic Use