

Product Name: Recombinant Human POLR2A Protein
Catalog #: PEH3545



Summary

Name	RPB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P24928
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IAEKINAGFG DDLNCIFNDD NAEKLVLRIR IMNSDENKMQ EEEEVVDKMD DDVFLRCIES NMLTDMTLQG IEQISKVYMH LPQTDNKKKI IITEDGEFKA LQEWILETDG VSLMRVLSEK

image



Background

Alternative Names	DNA-directed RNA polymerase II subunit RPB1 (RNA polymerase II subunit B1) (EC 2.7.7.6) (DNA-directed RNA polymerase II subunit A) (DNA-directed RNA polymerase III largest subunit) (RNA-directed RNA polymerase II subunit RPB1) (EC 2.7.7.48)
--------------------------	--

Note

Product Name: Recombinant Human POLR2A Protein
Catalog #: PEH3545



For Research Use Only , Not for Diagnostic Use