


Product Name: Recombinant Human LAMA1 Protein
Catalog #: PEH3546



Summary

Name	LAMA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P25391
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IERLQMSITE IMEKTTLNQT LDEDFLLPNS TLQNMQQNGT SLLEIMQIRD FTQLHQNATL ELKAAEDLLS QIQENYQKPL EELEVLKEAA SHVLSKHNNE LKAAEALVRE AEAKMQESNH LLLMVNANLR EFSDKKLHVQ EEQNLTSELI VQGRGLIDAA

image



Background

Alternative Names	Laminin subunit alpha-1 (Laminin A chain) (Laminin-1 subunit alpha) (Laminin-3 subunit alpha) (S-laminin subunit alpha) (S-LAM alpha)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LAMA1 Protein
Catalog #: PEH3546

