


Product Name: Recombinant Human ARNT Protein
Catalog #: PEH3550



Summary

Name	ARNT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P27540
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRQQQQQTELD M V P G R D G L A S Y N H S Q V V Q P V T T T G P E H S K P L E K S D G L F A Q D R D P R F S E I Y H N I N A D Q S K G I S S T V P A T Q Q L F S Q G N T F P P T P R P A E N F

image



Background

Alternative Names	Aryl hydrocarbon receptor nuclear translocator (ARNT protein) (Class E basic helix-loop-helix protein 2) (bHLHe2) (Dioxin receptor, nuclear translocator) (Hypoxia-inducible factor 1-beta) (HIF-1-beta) (HIF1-beta)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use