



Product Name: Recombinant Human PDIA3 Protein
Catalog #: PEH3554



Summary

Name	PDIA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30101
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEIPVVAIRT AKGEKFVMQE EFSRDGKALE RFLQDYFDGN LKRYLKSEPI PESNDGPVKV VVAENFDEIV NNENKDVLE FYAPWCGHCK NLEPKYKELG EKLSKDPNIV IAKMDATAND VSPYEVGRF PTIYFSPANK KLNPKKYEGG RELSDFISYL QREATNPPVI

image



Background

Alternative Names	Protein disulfide-isomerase A3 (EC 5.3.4.1) (58 kDa glucose-regulated protein) (58 kDa microsomal protein) (p58) (Disulfide isomerase ER-60) (Endoplasmic reticulum resident protein 57) (ER protein 57) (ERp57) (Endoplasmic reticulum resident protein 60) (ER protein 60) (ERp60)
--------------------------	--

Product Name: Recombinant Human PDIA3 Protein
Catalog #: PEH3554



Note

For Research Use Only , Not for Diagnostic Use