


Product Name: Recombinant Human LRPAP1 Protein
Catalog #: PEH3557



Summary

Name	AMRP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30533
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEEFRMEKLN QLWEKAQRLH LPPVRLAELH ADLKIQRDE LAWKKLKLDG LDEDGEKEAR LIRNLNVILA KYGLDGKKDA RQVTSNSLSG TQEDGLDDPR LEKLWHKAKT SGKFSGEELD KLRWREFLHHK EKVHEYNVLL ETLRTEEIH ENVISPSDLS DIKGSVLHSR HTELKEKLRV INQGLDRLRR VSHQGYSTEA EFEEPRVIDL WDLAQSANLT DKELEAFREE LKHFEAKIEK HNHYQKQLEI AHEKLRHAES VGDGERVSRV

image



Background

Alternative Names Alpha-2-macroglobulin receptor-associated protein (Alpha-2-MRAP) (Low density lipoprotein receptor-related protein-associated protein 1) (RAP)

Product Name: Recombinant Human LRPAP1 Protein
Catalog #: PEH3557



Note

For Research Use Only , Not for Diagnostic Use