


**Product Name: Recombinant Human PRKAR1B Protein**  
**Catalog #: PEH3558**



## Summary

<b>Name</b>	KAP1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P31321
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MASPPACPSE EDESLKGCEL YVQLHGIQQV LKDCIVHLCI SKPERPMKFL REHFEKLEKE ENRQILARQK SNSQSDSHDE EVSPTPPNPV VKARRRRGGV SAEVYTEEDA VSYVRKVIK DYKTMTALAK AISKNVLFAH LDDNERSDIF DAMFPVTHIA GETVIQQGNE GDNFYVVDQG EVDVYVNGEW VTNISEGGSF GELALIYGTP RAATVKAKTD LKLWGIDRDS YRRILMGSTL RKRKMYEFL SKVSILESLE KWERTLVADA LEPVQFEDGE KIVVQGE PGD DFYIITEGTA SVLQRRSPNE EYVEVGR LGP SDYFGEIALL LNRPRAATVV ARGPLKCVKL DRPRFERVLG PCSEILKRNI QRYNSFISLT V

## image



## Background

**Alternative Names** cAMP-dependent protein kinase type I-beta regulatory subunit

**Product Name: Recombinant Human PRKAR1B Protein**  
**Catalog #: PEH3558**

---



**Note**

For Research Use Only , Not for Diagnostic Use