


Product Name: Recombinant Human GBP1 Protein
Catalog #: PEH3559



Summary

Name	GBP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P32455
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VNAISSGDLP CMENAVLALA QIENSAAVQK AIAHYEQQMG QKVQLPTETL QELDLHRDS EREAIEVFIR SSFKDVDHLF QKELAAQLEK KRDDFCKQNNQ EASSDRCSAL LQVIFSPLEE EVKAGIYSKP GGYRLFVQKL QDLKKKYEE PRKGIQAEEI LQTYLKSKE MTDAILQTDQ TLTEKEKEIE VERVKAESAQ ASAKMLQEMQ RKNEQMMEQK ERSYQEHLKQ LTEKMENDRV QLLKEQERTL ALKLQEQEQL LKEGFQKESR

image



Background

Alternative Names	Guanylate-binding protein 1 (EC 3.6.5.-) (GTP-binding protein 1) (GBP-1) (HuGBP-1) (Guanine nucleotide-binding protein 1) (Interferon-induced guanylate-binding protein 1)
--------------------------	--

Product Name: Recombinant Human GBP1 Protein
Catalog #: PEH3559



Note

For Research Use Only , Not for Diagnostic Use