



Product Name: Recombinant Human TXLNA Protein
Catalog #: PEH3564



Summary

Name	TXLNA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P40222
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKQSNPKSSP GQPEAGPEGA QERPSQAAPA VEAEGPGSSQ APRKPEGAQA RTAQSGALRD VSEELSRQLE DILSTYCVDN NQGGPGEDGA QGEPAPEDA EKSRTYVARN GEPEPTPVVN GEKEPSKGDP NTEEIRQSDE VGDRDHRRPQ

image



Background

Alternative Names Alpha-taxilin

Note

For Research Use Only , Not for Diagnostic Use