


Product Name: Recombinant Human LIFR Protein
Catalog #: PEH3571



Summary

Name	LIFR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P42702
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CETHDLKEII CSWNPGRVTA LVGPRATSYT LVESFSGKYV RLKRAEAPTN ESYQLLFQML PNQEIYNFTL NAHNPLGRSQ STILVNITEK VYPHTPTSFK VKDINSTAVK LSWHLPGNFA KINFLCEIEI KKSNSVQEQR NVTIKGVENS SYLVALDKLN PYTLYTFRIR CSTETFWKWS KWSNKKQHLT TEASPSKGPD TWREWSSDGK NLIYWKPPLP INEANGKILS YNVSCSSDEE TQSLSEIPDP QHKAEIRLDK NDYIISVVAK NSVGSSPPSK

image



Background

Alternative Names Leukemia inhibitory factor receptor (LIF receptor) (LIF-R) (CD antigen CD118)

Product Name: Recombinant Human LIFR Protein
Catalog #: PEH3571



Note

For Research Use Only , Not for Diagnostic Use