

Product Name: Recombinant Human THOP1 Protein
Catalog #: PEH3577



Summary

Name	THOP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P52888
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MADAASPCSV VNDLRWDLSA QQIEERTREL IEQTKRVYDQ VGTQEFEDVS YESTLKALAD VEVTYTVQRN ILDFPQHVSP SKDIRTASTE ADKKLSEFDV EMSMREDVYQ RIVWLQEKVQ

image



Background

Alternative Names Thimet oligopeptidase (EC 3.4.24.15) (Endopeptidase 24.15) (MP78)

Note

For Research Use Only , Not for Diagnostic Use