


Product Name: Recombinant Human NCS1 Protein
Catalog #: PEH3581



Summary

Name	NCS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P62166
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVVEELTRKT YFTEKEVQQW YKGFIDKDCPS GQLDAAGFQK IYKQFFPFGD PTKFATFVFN VFDENKDGRI EFSEFIQALS VTSRGTLDK LRWAFKLYDL DNDGYITRNE MLDIVDAIQ MVGNTVELPE EENTPEKRVD RIFAMMDKNA DGKLTQEFQ EGSKADPSIV QALSLYDGLV

image



Background

Alternative Names	Neuronal calcium sensor 1 (NCS-1) (Frequenin homolog) (Frequenin-like protein) (Frequenin-like ubiquitous protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NCS1 Protein
Catalog #: PEH3581

