


Product Name: Recombinant Human ADAM17 Protein
Catalog #: PEH3589



Summary

Name	ADA17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78536
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QYESLSLFHP SNVEMLSSMD SASVRIIKPF PAPQTPGRLQ PAPVIPSAPA APKLDHQRMD TIQEDPSTDS HMDEDGFEEKD PFPNSSTAAK SFEDLTDHPV TRSEKAASFK LQRQNRVDSK

image



Background

Alternative Names	Disintegrin and metalloproteinase domain-containing protein 17 (ADAM 17) (EC 3.4.24.86) (Snake venom-like protease) (TNF-alpha convertase) (TNF-alpha-converting enzyme) (CD antigen CD156b)
--------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ADAM17 Protein
Catalog #: PEH3589

