


Product Name: Recombinant Human TOP2B Protein
Catalog #: PEH3594



Summary

Name	TOP2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02880
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LWSLTKEKVE ELIKQRDAKG REVNDLKRKS PSDLWKEDLA AFVEELDKVE SQEREDVLG MSGKAIKGKV GKPKVKKLQL EETMPSYGR RIPEITAMK ADASKLLKK KKGDLDTAAV KVEFDEEFSG

image



Background

Alternative Names DNA topoisomerase 2-beta (EC 5.6.2.2) (DNA topoisomerase II, beta isozyme)

Note

For Research Use Only , Not for Diagnostic Use