


Product Name: Recombinant Human PTH1R Protein
Catalog #: PEH3595



Summary

Name	PTH1R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03431
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VMTKEEQIFL LHRAQAQCEK RLKEVLQRPA SIMESDKGWT SASTSGKPRK DKASGKLYPE SEEDKEAPTG SRYRGRPCLP

image



Background

Alternative Names Parathyroid hormone/parathyroid hormone-related peptide receptor (PTH/PTHrP type I receptor) (PTH/PTHr receptor) (Parathyroid hormone 1 receptor) (PTH1 receptor)

Note

For Research Use Only , Not for Diagnostic Use