


Product Name: Recombinant Human ZBTB16 Protein
Catalog #: PEH3597



Summary

Name	ZBT16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q05516
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTLCDVVIMV DSQEFHAHRT VLACTSKMFE ILFHRNSQHY TLDLSPKTF QQILEYAYTA TLQAKAEDLD DLLYAAEILE IEYLEEQCLK MLETIQASDD NDTEATMADG GAEEEEEDRKA RYLKNIFISK

image



SDS



WB

Background

Alternative Names	Zinc finger and BTB domain-containing protein 16 (Promyelocytic leukemia zinc finger protein) (Zinc finger protein 145) (Zinc finger protein PLZF)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use