



Product Name: Recombinant Human EP300 Protein
Catalog #: PEH3609



Summary

Name	EP300_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q09472
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSNPPSTSST EVNSQAIAEK QPSQEVKMEA KMEVDQPEPA DTQPEDISES KVEDCKMEST ETEERSTELK TEIKEEDQP STSATQSSPA PGQSKKKIFK PEELRQALMP TLEALYRQDP ESLPFRQPVD PQLLGIPDYF DIVKSPMDLS TIKRKLDTGQ YQEPWQYVDD IWLMFNNAWL YNRKTSRVYK

image



Background

Alternative Names	Histone acetyltransferase p300 (p300 HAT) (EC 2.3.1.48) (E1A-associated protein p300) (Histone butyryltransferase p300) (EC 2.3.1.-) (Histone crotonyltransferase p300) (EC 2.3.1.-) (Protein 2-hydroxyisobutyryltransferase p300) (EC 2.3.1.-) (Protein propionyltransferase p300) (EC 2.3.1.-)
--------------------------	--

Product Name: Recombinant Human EP300 Protein
Catalog #: PEH3609



Note

For Research Use Only , Not for Diagnostic Use