


Product Name: Recombinant Human STX4 Protein
Catalog #: PEH3612



Summary

Name	STX4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12846
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERVALVVHPG TARLGSPDEE FFHKVRTIRQ TIVKLGNKVQ ELEKQQVTIL ATPLPEESMK QELQNLNRLDEI KQLGREIRLQ LKAIEPQKEE ADENYNYSVNT RMRKTQHGVL SQQFVELINK CNSMQSEYRE KNVERIRRQL KITNAGMVSD EELEQMLDSG QSEVFVSNIL KDTQVTRQAL NEISARHSEI QQLERSIREL HDIFTFLATE VEMQGEMINR IEKNILSSAD YVERGQEHVK

image



Background

Alternative Names Syntaxin-4 (Renal carcinoma antigen NY-REN-31)

Note

Product Name: Recombinant Human STX4 Protein
Catalog #: PEH3612



For Research Use Only , Not for Diagnostic Use