


Product Name: Recombinant Human GAB1 Protein
Catalog #: PEH3618



Summary

Name	GAB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13480
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QDYLLINCQ SKKPEPTRTH ADSAKSTSSE TDCNDNVPSH KNPASSQSKH GMNGFFQQQM IYDSPPSRAP SASVDSSLYN LPRSYSHDVL PKVSPSSTE A

image



Background

Alternative Names	GRB2-associated-binding protein 1 (GRB2-associated binder 1) (Growth factor receptor bound protein 2-associated protein 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use