


Product Name: Recombinant Human BIRC2 Protein
Catalog #: PEH3619



Summary

Name	BIRC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13490
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RYPHLLLEQLL STSDTTGEEN ADPPIHFGP GESSSEDAVM MNTPVVKSAL EMGFNRDLVK QTVQSKILTT GENYKTVNDI VSALLNAEDE KREEEKEKQA EEMASDDL SL IRKNRMALFQ QLTCVLPILD NLLKANVINK QEHDIIKQKT QIPLQARELI DTILVKGNA A ANIFKNCLKE IDSTLYKNLF VDKNMKYIPT EDVSGLSLEE QLRRRLQEERT

image



Background

Alternative Names	Baculoviral IAP repeat-containing protein 2 (EC 2.3.2.27) (Cellular inhibitor of apoptosis 1) (C-IAP1) (IAP homolog B) (Inhibitor of apoptosis protein 2) (hIAP-2) (hIAP2) (RING finger protein 48) (RING-type E3 ubiquitin transferase BIRC2) (TNFR2-TRAF-signaling complex protein 2)
--------------------------	---

Product Name: Recombinant Human BIRC2 Protein
Catalog #: PEH3619



Note

For Research Use Only , Not for Diagnostic Use