


Product Name: Recombinant Human MED21 Protein
Catalog #: PEH3620



Summary

Name	MED21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13503
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEEYAQLFAA LIARTAKDID VLIDSLPSEE STAALQAASL YKLEENHEA ATCLEDVVYR GDMILLEKIQS ALADIAQSQL

image



Background

Alternative Names	Mediator of RNA polymerase II transcription subunit 21 (Mediator complex subunit 21) (RNA polymerase II holoenzyme component SRB7) (RNAPII complex component SRB7) (hSrb7)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use