


Product Name: Recombinant Human PTK2B Protein
Catalog #: PEH3625



Summary

Name	FAK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14289
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEDFIQPSSR EEAQQLWEAE KVKMRQILDK QKQQMVEDYQ WLRQEEKSLD PMVYMNDKSP LTPEKEVGYL

image



Background

Alternative Names	Protein-tyrosine kinase 2-beta (EC 2.7.10.2) (Calcium-dependent tyrosine kinase) (CADTK) (Calcium-regulated non-receptor proline-rich tyrosine kinase) (Cell adhesion kinase beta) (CAK-beta) (CAKB) (Focal adhesion kinase 2) (FADK 2) (Proline-rich tyrosine kinase 2) (Related adhesion focal tyrosine kinase) (RAFTK)
--------------------------	---

Note

Product Name: Recombinant Human PTK2B Protein
Catalog #: PEH3625



For Research Use Only , Not for Diagnostic Use