


Product Name: Recombinant Human ZNF174 Protein
Catalog #: PEH3631



Summary

Name	ZN174_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15697
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PKQWVAVCMQ GQKVLEKTG SQLGEQELPD FQPQTPRRDL RESSPAEPSQ AGAYDRLSPH HWEKSPLLQE PTPKLAGTEA PRMRSDNKEN PQQEGAKGAK PCAVSAGRSK GNGLQNPEPR GANMSEPRLS RRQVSSPNAQ

image



Background

Alternative Names	Zinc finger protein 174 (AW-1) (Zinc finger and SCAN domain-containing protein 8)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use