



Product Name: Recombinant Human ENOX2 Protein
Catalog #: PEH3635



Summary

Name	ENOX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16206
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEEERLRPPS PPPVVHYSDH ECSIVAEKLK DDSKFSEAVQ TLLTWIERGE VNRRSANNFY SMIQSANS HV RRLVNEKAAH EKDMEEAKEK FKQALSGILI QFEQIVAVYH SASKQKAWDH FTKAQRKNIS VWCKQAEIR NIHNDELMGI RREEEMEMSD DEIEEMTETK ETEESALVSQ AEALKEENDS LRWQLDAYRN EVLLKQEQG KVRHREDDPNK EQQLKLLQQA LQGMQQHLLK VQEEYKKKEA ELEKLDDKL QVEKMLENLK EKESCASRLC ASNQDSEYPL EKTMNSSPIK

image



Background

Alternative Names	Ecto-NOX disulfide-thiol exchanger 2 (APK1 antigen) (Cytosolic ovarian carcinoma antigen 1) (Tumor-associated hydroquinone oxidase) (tNOX) [Includes: Hydroquinone [NADH] oxidase (EC 1.-.-.-); Protein disulfide-thiol
--------------------------	---

Product Name: Recombinant Human ENOX2 Protein
Catalog #: PEH3635



oxidoreductase (EC 1.-.-.-)]

Note

For Research Use Only , Not for Diagnostic Use