


Product Name: Recombinant Human CNGA2 Protein
Catalog #: PEH3636



Summary

Name	CNGA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16280
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MYIIKEGKLA VVADDGVTQY ALLSAGSCFG EISILNIKGS KMGNRRTANI RSLGYSDLFC LSKDDLMEAV TEYDPAKKVL EERGREILMK EGLLDENEVA TSMEVDVQEK LGQLETNMET LYTRFGRLLA EYTGAQQQKLK QRITVLETKM KQNNEDDYLSDGMNSPELAA

image



Background

Alternative Names	Cyclic nucleotide-gated olfactory channel (Cyclic nucleotide-gated cation channel 2) (Cyclic nucleotide-gated channel alpha-2) (CNG channel alpha-2) (CNG-2) (CNG2)
--------------------------	---

Note

Product Name: Recombinant Human CNGA2 Protein
Catalog #: PEH3636



For Research Use Only , Not for Diagnostic Use