


Product Name: Recombinant Human LAMA3 Protein
Catalog #: PEH3639



Summary

Name	LAMA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16787
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEKSQKEYEK LAASLNEARQ ELSDKVRELS RSAGKTSLVE EAEKHARSLQ ELAKQLEEIK RNASGDELVR CAVDAATAYE NILNAIKAAE DAANRAASAS ESALQTVIKE DLPRKAKTLS SNSDKLLNEA KMTQKCLKQE VSPALNNLQQ

image



Background

Alternative Names	Laminin subunit alpha-3 (Epiligrin 170 kDa subunit) (E170) (Epiligrin subunit alpha) (Kalinin subunit alpha) (Laminin-5 subunit alpha) (Laminin-6 subunit alpha) (Laminin-7 subunit alpha) (Nicein subunit alpha)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LAMA3 Protein
Catalog #: PEH3639

