


Product Name: Recombinant Human COL28A1 Protein
Catalog #: PEH3644



Summary

Name	COSA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2UY09
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKGEQGPQGF PGPKGTMGHG LPGQKGEHGE RGDVGGKGDK GEIGEPGSPG KQGLQGPKGD LGLTKEEIK LITEICGCGP KKETPLELV FVIDSSESVG PENFQIIKNF VKTMADRVAL DLATARIGII NYSHKVEKVA NLKQFSSKDD FKLAVDNMQY LGEPTYTATA LQAANDMFED ARPGVKKVAL VITDGQTDSR DKEKLTEVVK NASDTNVEIF VIGVVKKNDP NFEIFHKEMN LIATDPEHVY QFDDFFTLQD TLKQKLFQKI CEDFDSYLVQ IFGSSSPQPG FGMSGEELSE

image



Background

Alternative Names Collagen alpha-1(XXVIII) chain

Product Name: Recombinant Human COL28A1 Protein
Catalog #: PEH3644



Note

For Research Use Only , Not for Diagnostic Use