


Product Name: Recombinant Human CNST Protein
Catalog #: PEH3659



Summary

Name	CNST_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PJW8
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQIEPQDGCH PGDSVERSVT CLPSASDENE NQLDGDGHEH LTSSDSAMGK PQVSEQDSL N NESCTLSCE VAAGENLQNT LCEASRDEQA FLGKDKKIPG KRSPRSKKG T AKKIPPLGFS GDIAPLMQEK VLSAVTYAVD DEEAAEVNAN EQPEAPKLVL QSLFSLIRGE VEQLDSRALP LCLHQIAESY FQEEDYEKAM KFIQLERLYH EQLLANLSAI QEQWETKWKT VQPHTVTALR NSEKGFNGED FERLTKICAT HQDPLLSKHK IAAVEKSQER KCSTQLLVSE DPKEGGATTK ESESKTCLGT ESSKESQHTV EPLGSSPCCH QMDVQTDSPS LSVTAGKDHM

image



Background

Alternative Names Consortin

Product Name: Recombinant Human CNST Protein
Catalog #: PEH3659



Note

For Research Use Only , Not for Diagnostic Use