


Product Name: Recombinant Human SPESP1 Protein
Catalog #: PEH3660



Summary

Name	SPESP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UW49
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTTLDKSTGI GISTESDVP QLSGETAIEK PEEFGKHPES WNDDILKKI LDINSQVQQA LLSDTSNPAY REDIEASKDH LKRSLALAAA

image



Background

Alternative Names	Sperm equatorial segment protein 1 (SP-ESP) (Equatorial segment protein) (ESP) (Glycosylated 38 kDa sperm protein C-7/8)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use