


Summary

Name	AB1IP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z5R6
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RTIQAQKESL QNQHHSASLQ ASIFSGAASL GYGTVVAATG ISQYEDDLPP PPADPVLDLP LPPPPPEPLS QEEEEAQAKA DKIKLALAKL KEAKVKKLVV KVHMNDNSTK SLMVDERQLA RDVLDNLFK THCDCNVDWC LYEIYPELQI ERFFEDHENV VEVLSDWTRD TENKILFLEK

image



Background

Alternative Names	Amyloid beta A4 precursor protein-binding family B member 1-interacting protein (APBB1-interacting protein 1) (Proline-rich EVH1 ligand 1) (PREL-1) (Proline-rich protein 73) (Rap1-GTP-interacting adapter molecule) (RIAM) (Retinoic acid-responsive proline-rich protein 1) (RARP-1)
--------------------------	---

Product Name: Recombinant Human APBB1IP Protein
Catalog #: PEH3667



Note

For Research Use Only , Not for Diagnostic Use