


Product Name: Recombinant Human GALNT5 Protein
Catalog #: PEH3671



Summary

Name	GALT5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z7M9
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVKISVHMGR VSLKQEPKRS HSPSSDTSKL AAERDLNVTI SLSTDRPKQR SQAVANERAH PASTAVPKSG EAMALNKTKT QSKEVNANKH KANTSLPFPK FTVNSNRLRK

image



Background

Alternative Names	Polypeptide N-acetylgalactosaminyltransferase 5 (EC 2.4.1.41) (Polypeptide GalNAc transferase 5) (GalNAc-T5) (pp-GaNTase 5) (Protein-UDP acetylgalactosaminyltransferase 5) (UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 5)
--------------------------	--

Note

Product Name: Recombinant Human GALNT5 Protein
Catalog #: PEH3671



For Research Use Only , Not for Diagnostic Use