


Product Name: Recombinant Human MTDH Protein
Catalog #: PEH3672



Summary

Name	LYRIC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UE4
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRRSPPRKRE EAAAVPAAAP DDLALLKNLR SEEQKKKNRK KLSEKPKPNG RTVEVAEGEA VRTPQSVTAK QPPEIDKKNE KSKKNKKKSK SDAKAVQNSS RHDGKEVDEG

image



Background

Alternative Names	Protein LYRIC (3D3/LYRIC) (Astrocyte elevated gene-1 protein) (AEG-1) (Lysine-rich CEACAM1 co-isolated protein) (Metadherin) (Metastasis adhesion protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MTDH Protein
Catalog #: PEH3672

