

Product Name: Recombinant Human SESTD1 Protein
Catalog #: PEH3673



Summary

Name	SESD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86VW0
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQHSEWFAVY VELNQQAAL LNAGDEEDLV ELKSLQQQLS DVCYRQASQL EFRQNLLQAA LEFHGVAQDL SQQLDGLLGM LCVDVAPADG ASIQQTLKLL EEKLKSVDVG LQGLREKGQG LLDQISNQAS WAYGKDVITIE NKENVDHIQG VMEDMQLRKQ RCEDMVDVRR LKMLQMVQLF KCEEDAAQAV EWLSELLDAL LKTHIRLGDD

image



Background

Alternative Names SEC14 domain and spectrin repeat-containing protein 1 (Huntingtin-interacting protein-like protein) (Protein Solo)

Note

Product Name: Recombinant Human SESTD1 Protein
Catalog #: PEH3673



For Research Use Only , Not for Diagnostic Use