


Product Name: Recombinant Human VPS37A Protein
Catalog #: PEH3686



Summary

Name	VP37A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NEZ2
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AS IPTSQNGF GYKMPDVPDA FPELSELSVS QL TDMNEQEE VLLEQFLTLP QLKQIITDKD DLVKSIEELA RKNLLEPSL EAKRQTVLDK YELLTQMKST FEKKMQRQHE LSESCSASAL QARLKVAAHE AEEESDNIAE DFLEGKMEID DFLSSFMEKR TICHCRRAKE EKLQQAIA MH

image



Background

Alternative Names Vacuolar protein sorting-associated protein 37A (hVps37A) (ESCRT-I complex subunit VPS37A) (Hepatocellular carcinoma-related protein 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VPS37A Protein
Catalog #: PEH3686

