

Product Name: Recombinant Human LILRB4 Protein
Catalog #: PEH3688



Summary

Name	LIRB4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NHJ6
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAQRQADFQR PPGAAEPEPK DGGLQRRSSP AADVQGENFC AAVKNTQPED GVEMDTRQSP HDEDPQAVTY AKVKHSRPRR EMASPPSPLS GEFLDTKDRQ AEEDRQMDTE AAASEAPQDV TYAQLHSFTL RQKATEPPPS QEGASPAEPS

image



Background

Alternative Names	Leukocyte immunoglobulin-like receptor subfamily B member 4 (CD85 antigen-like family member K) (Immunoglobulin-like transcript 3) (ILT-3) (Leukocyte immunoglobulin-like receptor 5) (LIR-5) (Monocyte inhibitory receptor HM18) (CD antigen CD85k)
--------------------------	--

Note

Product Name: Recombinant Human LILRB4 Protein
Catalog #: PEH3688



For Research Use Only , Not for Diagnostic Use