


Product Name: Recombinant Human LILRB4 Protein
Catalog #: PEH3689



Summary

Name	LIRB4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8NHJ6
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MQAGPLPKPT LWAEPGSVIS WGNSTIWCQ GTLEAREYRL DKEESPAPWD RQNPLEPKNK ARFSIPSMTE DYAGRYRCYY RSPVGWSQPS DPLELVMTGA YSKPTLSALP SPLVTSGKSV TLLCQSRSPM DTFLLIKERA AHPLLHLRSE HGAQQHQAEF PMSPVTSVHG GTYRCFSSHG FSHYLLSHPS DPLELIVSGS LEDPRPSPTR SVSTAAGPED QPLMPTGSVP

image



Background

Alternative Names	Leukocyte immunoglobulin-like receptor subfamily B member 4 (CD85 antigen-like family member K) (Immunoglobulin-like transcript 3) (ILT-3) (Leukocyte immunoglobulin-like receptor 5) (LIR-5) (Monocyte inhibitory receptor HM18) (CD antigen CD85k)
--------------------------	--

Product Name: Recombinant Human LILRB4 Protein
Catalog #: PEH3689



Note

For Research Use Only , Not for Diagnostic Use