


Product Name: Recombinant Human GADD45GIP1 Protein
Catalog #: PEH3691



Summary

Name	G45IP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TAE8
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQKRREREQH IAECMAKMPQ MIVNWQQQQR ENWEKAQADK ERRARLQAEA QELLYQVDP RSARFQELLQ DLEKKERKRL

image



Background

Alternative Names	Growth arrest and DNA damage-inducible proteins-interacting protein 1 (39S ribosomal protein L59, mitochondrial) (MRP-L59) (CKII beta-associating protein) (CR6-interacting factor 1) (CRIF1) (Mitochondrial large ribosomal subunit protein mL64) (Papillomavirus L2-interacting nuclear protein 1) (PLINP) (PLINP-1) (p53-responsive gene 6 protein)
--------------------------	--

Note

Product Name: Recombinant Human GADD45GIP1 Protein
Catalog #: PEH3691



For Research Use Only , Not for Diagnostic Use