


Product Name: Recombinant Human CANT1 Protein
Catalog #: PEH3700



Summary

Name	CANT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WVQ1
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ADLDTESRAQ EENTWFSYLK KGYLTLSDSG DKVAVEWDKD HGVLESHLAE KGRGMELSDL IVFNGKLYSV DDRTGVVYQI EGSKAVPWVI LSDGDGTVEK GFKAEWLAVK DERLYVGGLG KEWTTTTGDV VNENPEWVKV VGYKGSVDHE

image



Background

Alternative Names	Soluble calcium-activated nucleotidase 1 (SCAN-1) (EC 3.6.1.6) (Apyrase homolog) (Putative MAPK-activating protein PM09) (Putative NF-kappa-B-activating protein 107)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CANT1 Protein
Catalog #: PEH3700

