


Product Name: Recombinant Human WASF1 Protein
Catalog #: PEH3706



Summary

Name	WASF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92558
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DKRKEKRKQK QKNLDRPHEP EKVPRAPHDR RREWQKLAQG PELAEEDANL LHKHIEVANG PASHFETRPQ

image



Background

Alternative Names	Wiskott-Aldrich syndrome protein family member 1 (WASP family protein member 1) (Protein WAVE-1) (Verprolin homology domain-containing protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use