

Summary

Name	FIG4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92562
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ICAVNLKCLI VKKFHKYEEE IDIHNEFFRP YELSSFDDTF CLAMTSSARD FMPKTVGIDP SPFTVRKPDE TGKSVLGNKS NREEAVLQRK TAASAPPPPS EEAVSSSEED DSGTDREEEG SVSQRSTPVK MTDAGDSAKV TENVVQPMKE LYGINLSDGL SEEDFSIYSR FVQLGQSQHK QDKNSQQPCS RCDSDGVIKLT PISAFSQDNI YEVQPPRVDR

image



Background

Alternative Names	Polyphosphoinositide phosphatase (EC 3.1.3.-) (Phosphatidylinositol 3,5-bisphosphate 5-phosphatase) (SAC domain-containing protein 3)
--------------------------	---

Note

Product Name: Recombinant Human FIG4 Protein
Catalog #: PEH3707



For Research Use Only , Not for Diagnostic Use