


Product Name: Recombinant Human VT11A Protein
Catalog #: PEH3711



Summary

Name	VT11A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96AJ9
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFAVLTAET SKIARVRLP PDEKKQMVAN VEKQLEEAKE LLEQMDLEVR EIPPQSRGMY SNRMRSYKQE MGKLETDFKR SRIAYSDEVN NELLGDDGNS SENQRAHLLD NTERLERSR RLEAGYQIAV ETEQIGQEML ENLSHDREKI QRARERLRET DANLGKSSRI

image



Background

Alternative Names	Vesicle transport through interaction with t-SNAREs homolog 1A (Vesicle transport v-SNARE protein Vti1-like 2) (Vti1-rp2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VT11A Protein
Catalog #: PEH3711

