


Product Name: Recombinant Human CDK5RAP3 Protein
Catalog #: PEH3714



Summary

Name	CK5P3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JB5
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQVLPMLRFV QKRGNSTVYE WRTGTEPSVV ERPHLEELPE QVAEDAIDWG DFGVEAVSEG TDSGISAEAA GIDWGIFPES DSKDPGGDGI DWGDDAVALQ ITVLEAGTQA PEGVARGPDA LTLLEYTETR NQFLDELMEL EIFLAQRAVE LSEEADVLSV SQFQLAPAIL QGQTKEKMVT MVSVLEDLIG KLTSLLQLHL FMILASPRYV DRVTEFLQQK LKQSLLALK KELMVQKQQE ALEEQAALP KLDLLLEKTK ELQKLEADI SKRYSGRPVN

image



Background

Alternative Names	CDK5 regulatory subunit-associated protein 3 (CDK5 activator-binding protein C53) (LXXLL/leucine-zipper-containing ARF-binding protein) (Protein HSF-27)
--------------------------	--



Product Name: Recombinant Human CDK5RAP3 Protein
Catalog #: PEH3714

Note

For Research Use Only , Not for Diagnostic Use