


Product Name: Recombinant Human PUS7 Protein
Catalog #: PEH3717



Summary

Name	PUS7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PZ0
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKRGALVVED NDSGVPVEET KKQKLSECSL TKGQDGLQND FLSISEDVPR PPDTVSTGKG GKNSEAQLED EEEEEEDGLS EEEEEEESES FADMMKHGLT EADVGITKFV SSHQGFGSIL KERYSDFFVH EIGKDGRISH LNDLSIPVDE EDPSEDIFTV LTAEEKQRLE ELQLFKNKET SVAIEVIEDT KEKRTIIHQAIKSLFPGLET KTEDREGKKY IVAYHAAGKK ALANPRKHSW PKSRGSYCHF VLYKENKDTM DAINVLSKYL RVKPNIFSVM GTKDKRAITV

image



Background

Alternative Names Pseudouridylate synthase 7 homolog (EC 5.4.99.-)

Product Name: Recombinant Human PUS7 Protein
Catalog #: PEH3717



Note

For Research Use Only , Not for Diagnostic Use