


Product Name: Recombinant Human MPHOSPH8 Protein
Catalog #: PEH3719



Summary

Name	MPP8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99549
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSAQTPKGRR LSGEERGLWS TDSAEEDKET KRNESKEYQ KRHSDSKEEK GRKEPKGLKT LKEIRNAFDL FKLTPEEKND VSENNRKREE IPLDFKTIDD HKTKENKQSL KERRNTRDET DTWAYIAAEG DQEVLDVCQ ADENS DGRQQ ILSLGMDLQL EWMKLEDFQK HLDGKDENFA ATDAIPSNVL RDAVKNGDYI

image



Background

Alternative Names M-phase phosphoprotein 8 (Two hybrid-associated protein 3 with RanBPM) (Twa3)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MPHOSPH8 Protein
Catalog #: PEH3719

