


Product Name: Recombinant Human UBAC1 Protein
Catalog #: PEH3723



Summary

Name	UBAC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSL1
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGASATDEEA RDELTEIFKK IRRKREFRAD ARAVISLMEM GFDEKEVIDA LRVNNNQNA ACEWLLGDRK PSPEELDKGI

image



Background

Alternative Names	Ubiquitin-associated domain-containing protein 1 (UBA domain-containing protein 1) (E3 ubiquitin-protein ligase subunit KPC2) (Glioblastoma cell differentiation-related protein 1) (Kip1 ubiquitination-promoting complex protein 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human UBAC1 Protein
Catalog #: PEH3723

