


Product Name: Recombinant Human TSEN34 Protein
Catalog #: PEH3724



Summary

Name	SEN34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSV6
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LALTSFKRQQ EESFQEQSAL AAEARETRRQ ELLEKITEGQ AAKKQKLEQA SGASSSQEAG SSQAAKEDET SDGQASGEQE EAGPSSSQAG PSNGVAPLPR SALLVQLATA

image



Background

Alternative Names	tRNA-splicing endonuclease subunit Sen34 (EC 4.6.1.16) (Leukocyte receptor cluster member 5) (tRNA-intron endonuclease Sen34) (HsSen34)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use