



Product Name: Recombinant Human RSPO3 Protein
Catalog #: PEH3725



Summary

Name	RSPO3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXY4
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEWNPWSPCT KKGKTCGFKR GTETRVREII QHPSAKGNLC PPTNETRKCT VQRKKCQKGE RGKKGRERKR KKPNGGESKE AIPDSKSLES SKEIPEQREN KQQQKKRKVQ DKQKSVSVST

image



Background

Alternative Names	R-spondin-3 (Protein with TSP type-1 repeat) (hPWTSR) (Roof plate-specific spondin-3) (hRspo3) (Thrombospondin type-1 domain-containing protein 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use